



Prescribed for Progress™



Inclusive Fitness PRO2® Total Body Exerciser Product Specifications

Order Code:
PRO234-INT(USA)

FEATURES

| | |
|--|---|
| Total body solution | Use with total body, upper body or lower body for versatile exercise |
| Dependent motion | Allows stronger extremities to assist weaker ones |
| Knee-to-elbow motion | Provides optimized torso rotation |
| Bi-directional movement and resistance | Pedaling in forward and reverse provides versatility, helps prolong exercise and promotes reciprocal muscle balance |
| Low starting resistance | 6 watts |
| Resistance range | 191 levels of resistance - 20 levels adjustable in .1 increments |
| Resistance system | 3 phase combination generator & eddy current brake - largest watt range in the industry |
| Easy to transport | Integrated transport wheels |
| Heart rate monitoring | Telemetric (chest strap required) |

ACCESSIBILITY & ADJUSTABILITY

| | |
|---|--|
| Color-coded for clear visual cues | Adjustment points and moving parts are color-coded |
| Adjustable arm and leg cranks | 3 length settings |
| Step-through seating | Adjustable step-through allows easy entry and safe and easy access for all users |
| Direct wheelchair access | Seat slides off for direct wheelchair access - ADA compliant |
| Heavy-duty wheelchair ramps | For manual and powered wheelchair access. Provides increased stabilization |
| Low support boots | Keep the feet on the pedals with sturdy straps around the back, top, and front of foot |
| Premium seat | Seat adjusts up/down, forward/back, swivels and reclines |
| 360° Therapist Assist Pedals | Adjust the seat forward and back from any point around the machine |
| Seat Recline | Recline to any point within a 23 degree range |
| Seat Swivel | 360 degree swivel that locks at each 90 degrees for safety |
| Seat height adjustment | Select seat height anywhere from 18.5"-23" |
| Seat removal handle | Integrated handle on seat back makes it easy to remove the seat |
| Accessible adjustment levers | Reach all adjustment levers from both sides of the seat |
| Zero-sloping forward/back seat adjustment | Adjust forward and back on level monorail for optimum control and safe positioning |

INTELLI-FIT™ CONSOLE

| | |
|--|---|
| 7" LCD full color touch screen display | Touch screen with clear, simple screens and large, readable font. |
| Audible beep | Yes |
| Tactile buttons | Tactile buttons in addition to touch screen |
| Resistance levels | 191 levels of resistance - 20 levels adjustable in .1 increments |
| Display feedback | Time, level, distance, METs, heart rate, calories, RPM, watts |
| One touch quick start | Yes |
| Control METs - Constant Work Program | Set constant METs level and the workload remains constant throughout the program. |
| Iso-Strength Program | Isokinetic program for safe, accommodating resistance |
| Power Burst | 8-second burst of high intensity training |

INTELLI-FIT™ CONSOLE (cont.)

| | | | | |
|---|--|-----------------------|-----------|---------------------------|
| Cool down | Automatic cool down begins after each exercise is completed | | | |
| Administrative mode | User defined settings, metric or US units of measure and other custom settings | | | |
| User-defined console languages | 9 languages - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch | | | |
| Download results to USB flash drive | Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file | | | |
| CSAFE port and 8-volt DC power port | Yes | | | |
| USB charger | For phones and other electronic devices | | | |
| Water bottle holder, accessory tray, reading rack | Integrated | | | |
| Personal Cooling Fan | Integrated dual fans with three speeds | | | |
| Programs | Quick Start | Constant Work - METs | Hills | Heart Fit Test & Training |
| | Manual | Constant Work - Watts | Fit-Quik® | Power Fit Test & Training |
| | Heart Rate | Iso-Strength | Random | Stress Test |

TECHNICAL SPECIFICATIONS

| | |
|------------------------------|---|
| Power Requirements | Self-generating; low-voltage AC adapter, included |
| Maximum user weight capacity | 450 lbs. |
| Unit weight | 260 lbs. |
| Unit size | 61" L x 30" W x 62" H |
| Color | Frame: Charcoal Cover: Cool Gray |
| Manufacturing | Assembled in USA |

WARRANTY

| | |
|-------|---------|
| Parts | 3 years |
| Labor | 1 year |

ENHANCEMENTS

| | |
|-----------------------|--|
| Fit-Key® software | Optional Fit-Key® software makes it easy to create specialized programs and track data |
| Entertainment options | Yes |

ACCESSORIES

| | |
|-----------------------------|---|
| Leg Stabilizers | Provides comfort and stabilization for users with leg weakness or spasticity. Durable steel construction with padded calf cuff. Large foot bed with straps for secure foot placement. |
| Assist Gloves | Designed for individuals who require help gripping the hand cranks. Made of heavy-duty closed cell foam that resists moisture. Flexible, comfortable and easy to clean. |
| Straight Grips | For those who have spasticity, weakness, or tremors in their upper extremities. |
| Heavy-Duty Wheelchair Ramps | For manual and powered wheelchair access. Provides increased stabilization |
| Polar Chest Strap | Transmits heart rate directly to console - works on all SCIFIT products |
| Seat Belt | Adjustable seat belt for added security and stability. Easy to attach and remove. Storage pouch included. Fits Premium seat only. |

SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, PRO2, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

PRO2 Specs - PRO234-INT-USA - 2015 - v1