



Prescribed for Progress™

AC5000M

FORWARD & REVERSE TREADMILL

The AC5000M is ideal for rehab environments as it has a reverse speed function up to 6kph and a -3% decline for downhill training.

AC5000M with Cool Grey Covers Handrails



- Low starting speed of 0.5kph, adjustable in 0.1 increments
- Forward & reverse speed range from 0.5kph to 19kph (reverse up to 6kph)
- Elevation range from -3% to 12%
- Maintenance-free, shock absorbing deck
- Generous walking surface at 56cm x 157cm
- Side handrail switches for speed and elevation
- Dual emergency stop
- Optional treadmill step

SPECIFICATIONS

Specification	AC5000M
Platform & Belt	Extra wide with large low static belt (56cm x 157cm)
Deck	Eucalyptus - maintenance-free no lubrication
Side Switches	Yes
Speed Increments	0.5kph
Speed Range	0.5-19kph
Reverse Speed	up to 6kph
Elevation Range	-3% to 12%
Medical CE IIa Certification	No
Full Side Handrails	Standard
Treadmill Step	Optional
Programs	Manual • Heart Rate • Programme Profiles
Feedback	Speed • Elevation • Heart Rate • Time • Calories • Distance • Pace • METS
Power	3HP Motor with inverter, mains power required
Length	207cm
Width	80cm
Height	160cm
Weight	220kg
Max. User Capacity	250kg
Color	Silver powder coated frame and cool grey covers
Warranty	3 years parts

Optional Accessories

- Treadmill Step