

SCIFIT

DC1000 Treadmill

Scientific Solutions

- Short step up height at only 7".
- Low starting speed at .4 and ranges up to 10 mph with speed adjustable in just .1 increments.
- Elevation ranges from 0%-10% adjustable in .5 increments.
- Maintenance-free, shock absorbing, reversible deck.
- 3 HP DC motor.
- Meets current leakage protection requirements.
- Walking surface of 19 1/2" x 58 1/2".
- Emergency stop with magnetic clip and button on the console.
- Easy grip, extended handrails included.
- Front rollers for easy moving and transport.
- Contact heart rate located on front handlebar.
- Heart Rate compatible with most wireless technology.

Intelli-Fit[™] Features and Programs

- 7" LCD display with user friendly interface using tactile keypad, color coded with audible beep.
- Programs: One button Quick Start, Manual, Heart Rate (optional chest strap must be used), 7 pre-programmed elevation profiles, 7 pre-programmed speed profiles.
- Administrative mode: includes 6 user defined languages, metric or US units of measure and other custom settings.
- Built-in personal twin Fans: 3 levels of intensity at low, medium and high.
- Automatic cool down after each exercise is completed.
- Download results to USB flash drive. Simply plug your USB key in the Fit-Key port and your work-out will be copied in a text file for saving and printing.
- Display Feedback: Heart rate, time, speed, elevation, calories, distance in miles, pace, ascent (ft.) and METs.
- CSAFE port and 8-volt DC power port for TV or other controllers.
- USB charger for phones, mp3 players and other electronic devices.

Warranty

Limited warranty of only 6 hours of use per day: 3 year parts, 1 year labor with a 5 year frame and drive system warranty. Exceptions: 1 year warranty on: deck, running belt, contact heart rate grips, heart rate receiver and transmitter. Wear items (grips and rubber foot pads) are excluded

Optional Accessories



USB Fit-Keys



Heart Rate Chest Strap



Technical Specifications

- Overall Dimensions: 73"L x 32" W x 58" H
- Weight of unit: 222 lbs pre packaged
- Maximum user weight: 375 lbs
- 110v standard
- Electrical requirements: dedicated 15 amp circuit, grounded.

Order Code

DC1000-INT