



Scientific Solutions For Fitness™

Product Specifications

SCIFIT AC5000M Medical Treadmill

Features

- Speed range from .1 - 10 MPH with reverse speeds to 4 MPH.
- -3% to +12% elevation in .5 increments of adjustment.
- Side handrail switches for speed and elevation.
- Reversible eucalyptus deck with shock absorption.
- Commercial grade 3" crowned rollers.
- Oversized walking surface.
- Fit-Key™ technology, optional.
- Heart Rate compatible with most wireless technology, and contact heart rate.
- Integrated reading rack and water bottle holder.
- Dual emergency stop for added safety.
- Includes side hand rails.

Performance

- Powerful 3 HP AC motor with inverter.
- Speed adjustable in .1 MPH increments.

Programs

- Manual - One Button Quick Start.
- Heart Rate - Treadmill adjusts elevation to maintain target heart rate.
- Profile Programs incorporating changes in speed and elevation: Peak, Twin Peaks, Multi-Peak, Ramp, Climb Course, Progressive, and Speed Sprint.

Feedback

speed, elevation, distance, pace, time, calories, heart rate, METS

Dimensions

- Height 64"
- Length 81"
- Width 36.5" (with hand rails)
- Belt Size 62" x 22" with SCIFIT logo on belt for added safety
- Weight 520 lbs.
- User capacity 500 lbs.

Color

Frame - Silver powder coat with black accents
Covers - Black & Silver Metallic with UV Protection

Power

110V or 220V

Warranty

3 years parts, 1 year labor, exclude wear items

Order Code

AC5000M