



Prescribed for Progress™



DC1000 Treadmill Product Specifications

Order Code: DC1000-INT

FEATURES

Low starting speed	0.4 mph
Speed range	0.4 mph to 10 mph adjustable in .1 increments
Elevation range	0% to 10% adjustable in .5 increments
Low step-up height	7 inches
Multiple emergency stop	Safety clip, stop button by the console, and stop button on the touch screen provided added safety
Extended handrails, included	Provides extra support, safety and comfort for entering, exiting and during use
Easy to transport	Integrated transport wheels
Heart rate monitoring	Contact and telemetric (chest strap required)

BELT, DECK AND DRIVE SYSTEM

Walking surface	19.5" x 58.5"
Continuous 2-ply carbon fiber weave belt	Helps dissipate static electricity
Polyester synthetic backing	Provides a long lasting surface
Belt splices	Thermo form welded at specific angles to run smooth and quiet
High tensile strength rating	Overlapping finger splice design to provide the ultimate in tensile/tear strength
Eucalyptus hardwood deck	Strong and durable
Very low maintenance	No waxing or lubrication is required
Reversible deck	Extends the life of the deck
3HP DC motor	Provides smooth, natural feel

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Yes
Easy operation	Touch screen with clear, simple screens and large, readable font.
Audible beep	Yes
Tactile buttons	Tactile buttons in addition to touch screen
Display feedback	Time, speed, elevation, distance, ascent, pace, calories, METs, heart rate
One touch quick start	Yes
Cool down	Automatic cool down begins after each exercise is completed
Administrative mode	User defined settings, metric or US units of measure and other custom settings
Console Languages	9 languages available
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file
CSAFE port and 8-volt DC power port	Yes
USB charger	For phones and other electronic devices
Water bottle holder, accessory tray, reading rack	Integrated
Fan	Integrated dual fans with three speeds
Programs	Quick Start 7 Speed Profiles Manual 7 Elevation Profiles Heart Rate

TECHNICAL SPECIFICATIONS

Current Leakage Protection Standards	Conforms to safety requirements of Standard IEC60601-1:2005 and NFPA-99	
Certifications	UL1647; CAN/CSA-C22.2 No 68-92; EN60335	
Power Requirements	110v; dedicated 15 amp circuit, grounded	
Maximum user weight capacity	375 lbs.	
Unit weight	222 lbs.	
Unit size	73"L x 32"W x 58"H	
Color	Frame: Silver Cover: Black	
Manufacturing	Assembled in USA	

WARRANTY

Limited Warranty	Only 6 hours of use per day	
Parts: 3 years	Labor: 1 year	Frame & Drive System: 5 years

ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data	
Entertainment options	Yes	

ACCESSORIES

Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products	
-------------------	---	--



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

DC1000 Specs - DC1000-INT - 2015 - v1