



Prescribed for Progress™



Inclusive Fitness PRO2® Total Body Exerciser Product Specifications

Order Code:
PRO234-INT(USA)

FEATURES

Total body solution	Use with total body, upper body or lower body for versatile exercise
Dependent motion	Allows stronger extremities to assist weaker ones
Knee-to-elbow motion	Provides optimized torso rotation
Bi-directional movement and resistance	Pedaling in forward and reverse provides versatility, helps prolong exercise and promotes reciprocal muscle balance
Low starting resistance	6 watts
Resistance range	191 levels of resistance - 20 levels adjustable in .1 increments
Resistance system	3 phase combination generator & eddy current brake - largest watt range in the industry
Easy to transport	Integrated transport wheels
Heart rate monitoring	Telemetric (chest strap required)

ACCESSIBILITY & ADJUSTABILITY

Color-coded for clear visual cues	Adjustment points and moving parts are color-coded
Adjustable arm and leg cranks	3 length settings
Step-through seating	Adjustable step-through allows easy entry and safe and easy access for all users
Direct wheelchair access	Seat slides off for direct wheelchair access - ADA compliant
Heavy-duty wheelchair ramps	For manual and powered wheelchair access. Provides increased stabilization
Low support boots	Keep the feet on the pedals with sturdy straps around the back, top, and front of foot
Premium seat	Seat adjusts up/down, forward/back, swivels and reclines
360° Therapist Assist Pedals	Adjust the seat forward and back from any point around the machine
Seat Recline	Recline to any point within a 23 degree range
Seat Swivel	360 degree swivel that locks at each 90 degrees for safety
Seat height adjustment	Select seat height anywhere from 18.5"-23"
Seat removal handle	Integrated handle on seat back makes it easy to remove the seat
Accessible adjustment levers	Reach all adjustment levers from both sides of the seat
Zero-sloping forward/back seat adjustment	Adjust forward and back on level monorail for optimum control and safe positioning

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Touch screen with clear, simple screens and large, readable font.
Audible beep	Yes
Tactile buttons	Tactile buttons in addition to touch screen
Resistance levels	191 levels of resistance - 20 levels adjustable in .1 increments
Display feedback	Time, level, distance, METs, heart rate, calories, RPM, watts
One touch quick start	Yes
Control METs - Constant Work Program	Set constant METs level and the workload remains constant throughout the program.
Iso-Strength Program	Isokinetic program for safe, accommodating resistance
Power Burst	8-second burst of high intensity training

INTELLI-FIT™ CONSOLE (cont.)

Cool down	Automatic cool down begins after each exercise is completed			
Administrative mode	User defined settings, metric or US units of measure and other custom settings			
User-defined console languages	9 languages - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch			
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file			
CSAFE port and 8-volt DC power port	Yes			
USB charger	For phones and other electronic devices			
Water bottle holder, accessory tray, reading rack	Integrated			
Personal Cooling Fan	Integrated dual fans with three speeds			
Programs	Quick Start	Constant Work - METs	Hills	Heart Fit Test & Training
	Manual	Constant Work - Watts	Fit-Quik®	Power Fit Test & Training
	Heart Rate	Iso-Strength	Random	Stress Test

TECHNICAL SPECIFICATIONS

Power Requirements	Self-generating; low-voltage AC adapter, included
Maximum user weight capacity	450 lbs.
Unit weight	260 lbs.
Unit size	61" L x 30" W x 62" H
Color	Frame: Charcoal Cover: Cool Gray
Manufacturing	Assembled in USA

WARRANTY

Parts	3 years
Labor	1 year

ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data
Entertainment options	Yes

ACCESSORIES

Leg Stabilizers	Provides comfort and stabilization for users with leg weakness or spasticity. Durable steel construction with padded calf cuff. Large foot bed with straps for secure foot placement.
Assist Gloves	Designed for individuals who require help gripping the hand cranks. Made of heavy-duty closed cell foam that resists moisture. Flexible, comfortable and easy to clean.
Straight Grips	For those who have spasticity, weakness, or tremors in their upper extremities.
Heavy-Duty Wheelchair Ramps	For manual and powered wheelchair access. Provides increased stabilization
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products
Seat Belt	Adjustable seat belt for added security and stability. Easy to attach and remove. Storage pouch included. Fits Premium seat only.

SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, PRO2, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

PRO2 Specs - PRO234-INT-USA - 2015 - v1