



Prescribed for Progress™

ISO 1000R

FORWARD-ONLY RECUMBENT BIKE

With a true adjustable step through that adjusts up to a full 58cm, SCIFIT recumbent bikes offer direct wheelchair accessibility, along with, safe and easy access.

- Forward-only resistance
- Step-through easy entry on horizontal monorail for safe positioning
- Wheelchair Accessibility
- Custom 3-piece crank system with oversized pedals
- Iso-Strength safe, accommodating strength program
- Very low starting resistance
- Available with Premium or Standard Seat

ISO1000R Foward-only with Standard Seat and Charcoal Covers



SPECIFICATIONS

Specification	ISO1000R
Workout Position	Recumbent
Console	Intelli-Fit
Seat System	Wheelchair Accessibility • Standard Seat with forward/backward adjusters
Pedals	Oversized • adjustable footstraps
Cranks	Fixed
Medical CE IIa Certification	No
Bi-directional Resistance	No
Resistance	Workload levels range from 1 - 20, adjustable in 0.1 increments, providing 191 levels of resistance
Programs	Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training
Feedback	Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS
Power	Self-generating from just 6 watts with auto recharge and battery back up
Length	150cm
Width	66cm
Height	132cm
Weight	95kg
Max. User Capacity	227kg
Color	Silver powder coated frame and charcoal covers
Warranty	3 years parts

Seat Options

- Premium Seat (204kg max user weight)

Optional Accessories

- Adjustable Cranks
- Low Support Boots

