

ISO7000R

BI-DIRECTIONAL RECUMBENT BIKE

With the same true adjustable step through that adjusts up to a full 58cm with wheelchair accessibility, along with, safe and easy access as the ISO1000R, this model also offers bi-directional resistance.

- Bi-directional resistance for versatility and reciprocal muscle balance
- Step-through easy entry on horizontal monorail for safe positioning
- Wheelchair Accessibility
- Custom 3-piece crank system with oversized pedals
- Iso-Strength safe, accommodating strength program
- Very low starting resistance
- Available with Premium or Standard Seat

ISO7000R Bi-directional with Premium Seat and Cool Grey Covers



SPECIFICATIONS

| Specification | ISO7000R |
|------------------------------|--|
| Workout Position | Recumbent |
| Console | Intelli-Fit |
| Seat System | Wheelchair Accessibility • Premium Seat has 360-degree swivel, recline & full adjustability |
| Pedals | Oversized • adjustable footstraps |
| Cranks | Fixed |
| Medical CE IIa Certification | No |
| Bi-directional Resistance | Yes |
| Resistance | Workload levels range from 1 - 20, adjustable in 0.1 increments, providing 191 levels of resistance |
| Programs | Manual • Iso-Strength • Heart Rate • Constant Work • Random • Hill Profiles • Fit-Quik • Power Burst • Heart-Fit Test • Heart-Fit Training • Power-Fit Test • Power-Fit Training |
| Feedback | Heart Rate • Time • RPM • Calories • Distance • Level • Watts • METS |
| Power | Self-generating from just 6 watts with auto recharge and battery back up |
| Length | 150cm |
| Width | 66cm |
| Height | 132cm |
| Weight | 95kg |
| Max. User Capacity | 204kg |
| Color | Silver powder coated frame and cool grey covers |
| Warranty | 3 years parts |

Seat Options

- Standard Seat (227kg max user weight)

Optional Accessories

- Adjustable Cranks
- Low Support Boots