



Prescribed for Progress™



AC5000M Medical Treadmill Product Specifications

Order Code: AC5000M-INT

FEATURES

| | |
|-------------------------------|--|
| Low starting speed | 0.1 mph |
| Speed range | 0.1 mph to 12 mph adjustable in .1 increments |
| Reverse walking | Reversible speed goes up to 4 mph |
| Decline | Up to -3% decline for downhill training |
| Elevation range | -3% to 12% adjustable in .5 increments |
| Speed & elevation adjustments | Side handrail switches |
| Multiple emergency stop | Safety clip, stop button by the console, and stop button on the touch screen provided added safety |
| Extended handrails, included | Provides extra support, safety and comfort for entering, exiting and during use |
| Easy to transport | Integrated transport wheels |
| Heart rate monitoring | Contact and telemetric (chest strap required) |

BELT, DECK AND DRIVE SYSTEM

| | |
|--|---|
| Generous walking surface | 22" x 62" |
| Continuous 2-ply carbon fiber weave belt | Helps dissipate static electricity |
| Polyester synthetic belt backing | Provides a long lasting surface |
| Thermo form welded belt splices | Thermo form welded at specific angles to run smooth and quiet |
| High tensile strength rating belt | Overlapping finger splice design to provide the ultimate in tensile/tear strength |
| Eucalyptus hardwood deck | 1" deck is strong and durable |
| Maintenance-free deck | No waxing or lubrication is required |
| Reversible deck | Extends the life of the deck |
| Multiple zone cushioning system | 6 independent bumpers provide superior shock absorption |
| 3" crowned front and rear rollers | Helps maintain belt alignment and increase the life of the belt |
| 3HP AC motor | Maintenance-free motor with variable frequency torque vector drive system |
| Power Factor Correction System | Reduces power consumption by 25% |

INTELLI-FIT™ CONSOLE

| | |
|--|--|
| 7" LCD full color touch screen display | Yes |
| Easy operation | Touch screen with clear, simple screens and large, readable font. |
| Audible beep | Yes |
| Tactile buttons | Tactile buttons in addition to touch screen |
| Display feedback | Time, speed, elevation, distance, ascent, pace, calories, METs, heart rate |
| One touch quick start | Yes |
| Cool down | Automatic cool down begins after each exercise is completed |
| Administrative mode | User defined settings, metric or US units of measure and other custom settings |
| Console Languages | 9 languages available - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch |

INTELLI-FIT™ CONSOLE (cont.)

| | | |
|---|--|----------------------|
| Download results to USB flash drive | Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file | |
| CSAFE port and 8-volt DC power port | Yes | |
| USB charger | For phones and other electronic devices | |
| Water bottle holder, accessory tray, reading rack | Integrated | |
| Fan | Integrated dual fans with three speeds | |
| Programs | Quick Start | 7 Speed Profiles |
| | Manual | 7 Elevation Profiles |
| | Heart Rate | |

TECHNICAL SPECIFICATIONS

| | | |
|--------------------------------------|---|--|
| Current Leakage Protection Standards | Conforms to safety requirements of Standard IEC60601-1:2005 and NFPA-99 | |
| Certifications | UL1647; CAN/CSA-C22.2 No 68-92; EN60601; EN957-6; ISO 9001:2008; ISO 13485:2003/2012; NFPA-99 | |
| Power Requirements | 110v standard; dedicated 15 amp circuit, grounded (220v optional) | |
| Maximum user weight capacity | 550 lbs. | |
| Unit weight | 458 lbs. | |
| Unit size | 83" L x 34" W x 63" H | |
| Color | Frame: Charcoal Cover: Black & Cool Gray | |

WARRANTY

| | |
|----------------------|---------|
| Parts | 3 years |
| Labor | 1 year |
| Frame & Drive System | 5 years |

ENHANCEMENTS

| | |
|-----------------------|--|
| Fit-Key® software | Optional Fit-Key® software makes it easy to create specialized programs and track data |
| Entertainment options | Yes |

ACCESSORIES

| | |
|-------------------|--|
| Safety Step | 6" high step with non-slip safety strips and anti-tip design for easy entry and exit |
| Polar Chest Strap | Transmits heart rate directly to console - works on all SCIFIT products |



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

AC5000M Specs - AC5000M-INT - 2016 -v1