



Prescribed for Progress™



AC5000 Treadmill Product Specifications

Order Code: AC5000-INT

FEATURES

Low starting speed	0.1 mph
Speed range	0.1 mph to 12 mph adjustable in .1 increments
Elevation range	0% - 15% adjustable in .5 increments
Speed & elevation adjustments	Side handrail switches
Multiple emergency stop	Safety clip, stop button by the console, and stop button on the touch screen provided added safety
Extended handrails, optional	Provides extra support, safety and comfort for entering, exiting and during use
Easy to transport	Integrated transport wheels
Heart rate monitoring	Contact and telemetric (chest strap required)

BELT, DECK AND DRIVE SYSTEM

Generous walking surface	22" x 62"
Continuous 2-ply carbon fiber weave belt	Helps dissipate static electricity
Polyester synthetic belt backing	Provides a long lasting surface
Thermo form welded belt splices	Thermo form welded at specific angles to run smooth and quiet
High tensile strength rating belt	Overlapping finger splice design to provide the ultimate in tensile/tear strength
Eucalyptus hardwood deck	1" deck is strong and durable
Maintenance-free deck	No waxing or lubrication is required
Reversible deck	Extends the life of the deck
Multiple zone cushioning system	6 independent bumpers provide superior shock absorption
3" crowned front and rear rollers	Helps maintain belt alignment and increase the life of the belt
3HP AC motor	Maintenance-free motor with variable frequency torque vector drive system
Power Factor Correction System	Reduces power consumption by 25%

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Yes
Easy operation	Touch screen with clear, simple screens and large, readable font.
Audible beep	Yes
Tactile buttons	Tactile buttons in addition to touch screen
Display feedback	Time, speed, elevation, distance, ascent, pace, calories, METs, heart rate
One touch quick start	Yes
Cool down	Automatic cool down begins after each exercise is completed
Administrative mode	User defined settings, metric or US units of measure and other custom settings
Console Languages	9 languages available - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file
CSAFE port and 8-volt DC power port	Yes

INTELLI-FIT™ CONSOLE (cont.)

USB charger	For phones and other electronic devices	
Water bottle holder, accessory tray, reading rack	Integrated	
Fan	Integrated dual fans with three speeds	
Programs	Quick Start	7 Speed Profiles
	Manual	7 Elevation Profiles
	Heart Rate	

TECHNICAL SPECIFICATIONS

Current Leakage Protection Standards	Conforms to safety requirements of Standard IEC60601-1:2005 and NFPA-99	
Certifications	UL1647; CAN/CSA-C22.2 No 68-92; EN60601; EN957-6; ISO 9001:2008; ISO 13485:2003/2012; NFPA-99	
Power Requirements	110v standard; dedicated 15 amp circuit, grounded (220v optional)	
Maximum user weight capacity	550 lbs.	
Unit weight	458 lbs.	
Unit size	83" L x 34" W x 63" H	
Color	Frame: Charcoal	Cover: Black & Charcoal

WARRANTY

Parts	3 years
Labor	1 year
Frame & Drive System	5 years

ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data
Entertainment options	Yes

ACCESSORIES

Safety Step	6" high step with non-slip safety strips and anti-tip design for easy entry and exit
Extended Handrail Set	Provides extra support, safety and comfort for entering, exiting and during use
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

AC5000 specs - AC5000-INT - 2016 - v1