



Prescribed for Progress™



## ISO7000R Recumbent Bike Product Specifications

Order Code: ISO7027R-INT

### FEATURES

Bi-directional movement and resistance	Pedaling in forward and reverse provides versatility, helps prolong exercise and promotes reciprocal muscle balance
Low starting resistance	6 watts
Resistance range	191 levels of resistance - 20 levels adjustable in .1 increments
Resistance system	3 phase combination generator & eddy current brake - largest watt range in the industry
Large, self-righting pedals	Strap adjusts easily to fit the user
Handlebars	Handlebars on seat and console provide a variety of hand grip options
Easy to transport	Integrated transport wheels
Heart rate monitoring	Contact and telemetric (chest strap required)

### ACCESSIBILITY & ADJUSTABILITY



Step-through seating	Adjustable step-through allows easy entry and safe and easy access for all users
Direct wheelchair access	Seat slides off for direct wheelchair access - ADA compliant
Bariatric Seat	Wide seat adjusts forward and back, has side handles and 600 lb user weight capacity
360° Therapist Assist Pedals	Adjust the seat forward and back from any point around the bike
Seat removal handle	Integrated handle on seat back makes it easy to remove the seat
Zero-sloping forward/back seat adjustment	Adjust forward and back on level monorail for optimum control and safe positioning

### INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Touch screen with clear, simple screens and large, readable font.			
Audible beep	Yes			
Tactile buttons	Tactile buttons in addition to touch screen			
Resistance levels	191 levels of resistance - 20 levels adjustable in .1 increments			
Display feedback	Time, level, distance, METs, heart rate, calories, RPM, watts			
One touch quick start	Yes			
Control METs - Constant Work Program	Set constant METs level and the workload remains constant throughout the program.			
Iso-Strength Program	Isokinetic program for safe, accommodating resistance			
Power Burst	8-second burst of high intensity training			
Cool down	Automatic cool down begins after each exercise is completed			
Administrative mode	User defined settings, metric or US units of measure and other custom settings			
User-defined console languages	9 languages - English, Spanish, French, German, Russian, Japanese, Chinese, Italian, Dutch			
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file			
CSAFE port and 8-volt DC power port	Yes			
USB charger	For phones and other electronic devices			
Water bottle holder, accessory tray, reading rack	Integrated			
Personal Cooling Fan	Integrated dual fans with three speeds			
Programs	Quick Start	Constant Work - METs	Hills	Heart Fit Test & Training
	Manual	Constant Work - Watts	Fit-Quik®	Power Fit Test & Training
	Heart Rate	Iso-Strength	Random	Stress Test

## TECHNICAL SPECIFICATIONS

Power Requirements	Self-generating; low-voltage AC adapter, included
Maximum user weight capacity	600 lbs.
Unit weight	243 lbs.
Unit size	58" L x 24" W x 53" H
Color	Frame: Charcoal Cover: Cool Gray

## WARRANTY

Parts	3 years
Labor	1 year

## ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data
Entertainment options	Yes

## ACCESSORIES

Adjustable Pedal Cranks	Provides three adjustments (5", 6" and 7") to give the user versatility to adjust their range of motion
Low Support Boots	Keeps the foot on the pedal with sturdy straps around the back, top, and front of the foot.
Sports Performance Pedals	High-quality design for high-intensity training. Compatible with Shimano SPD cycling shoes. Reversible for cleat connection. Easily removable.
Wheelchair platform	Provides a sturdy foundation for wheelchair use
Heavy duty wheelchair ramps	For manual and powered wheelchair access. Provides increased stabilization
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time without notice and without incurring any obligation.

©Copyright 2015, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit and Prescribed for Progress are trademarks of SCIFIT Systems, Inc. Fit-Key Technology is Patented.

ISO7000R Specs - ISO7027R-INT - 2016 - v1